

Hospital/Institution Name: Radiological Practice Dres. Müller, Hildenbrand

City/Country: Köln Germany



Contact Persons (Phone, E-Mail): Herr Dr. T. Müller ☎ +49 221 - 312508

Consultant (Address, Phone, E-Mail):

Connected Modalities DICOM: 1 x CT, 1 x MRI, 1 x CR, 1 x NUC

Connected Modalities ANALOG:

Installation date: November 2003

Investment Volume ca.: 60.000 € - 80.000 €

Patients per day: GB per day: 50 1 GB

Server size: TB per year: 0,5 TB 0,3 TB

Number of Diagnosis Stations: 2

HIS/RIS Integration with Vendor: Datavital DavidX

Telecommunication with:

Special Applications or Integrations: Interface to Praxis management system. All modality with DICOM Worklist. DICOM Worklist come from VEPRO PACS

