



Hospital/Institution Name:	Kreiskrankenhaus Belzig GmbH		
City/Country:	Belzig	Germany	
Contact Persons (Phone, E-Mail):	CA Dr. Manfred Volkmann	+49 (0) 33841 - 93 - 100	
			✉ manfred.volkmann@kkh-belzig.de
Consultant (Address, Phone, E-Mail):	Hr. Torsten Grätz, CEO	+49 (0)333841 - 93 - 420	
			✉ torsten.graetz@kkh-belzig.de
Connected Modalities DICOM:	1 x CT, 1 x CR, 1x XA (Angio), 1 Laserimager		
Connected Modalities ANALOG:	-		
Installation date:	Feb. 2005		
Investment Volume ca.:	160.000 €		
Patients per day:	GB per day:	50	2 GB
Server size:	TB per year:	3 TB	0,6 TB
Number of Diagnosis Stations:	6		
HIS/RIS Integration with Vendor:	RIS (GAP/ iSOFT) with 10 workstation and HL7-Interface for report data		
Telecommunication with:	Klinikum "Ernst von Bergmann" – Postdam Home Office Dr. Seipt		
Special Applications or Integrations:	Clinic wide Image distribution via DICOM from the start, Clinic wide Report distribution via VEPRO PACS, HL7-Interface to RIS "iSOFT" Data import from a further PACS-solution All Modalities with DICOM Worklist		



Netzwerkplan KIS- RIS – PACS und Modalitäten



Außenansicht des KKH Belzig mit Hubschrauberlandeplatz über der Notaufnahme



Außenansicht Kreiskrankenhaus Belzig



CA Dr. Volkmann bei der Vorbereitung der Filmbesprechung



MTRA bei der Arbeit am CD-/DVD-Roboter



Mitarbeiter der DV-Abteilung bei Kontrolle am PACS-Server



logie beim Befunden an einer hochauflösenden Doppelmonitor-Station



inierter Arbeitsplatz für Filmbesprechung u. Notaufnahme



Auch an WEB- Arbeitsplätzen stehen alle Bilder als DICOM- Datensatz in Sekundenschnelle zur Verfügung