


Hospital/Institution Name:	Radiological Practice Dres. Reinheimer, Simon, Stölben		
City/Country:	Wittlich	Germany	
Contact Persons (Phone, E-Mail):	Hr. Dr. med. A. Stölben	☎ +49 6571 - 209 85	
	Hr. R. Rosenbaum	☎ +49 6571 - 151030	
Consultant (Address, Phone, E-Mail):			
Connected Modalities DICOM:	1 x CT, 1 x MRI, 1 x CR, 1 x DL, 1 x 3D Workstation		
Connected Modalities ANALOG:			
Installation date:	March 2004		
Investment Volume ca.:	400.000 € - 450.000 €		
Patients per day:	GB per day:	100	4 GB
Server size:	TB per year:	5,1 TB	0,5 TB
Number of Diagnosis Stations:	3 (each with 4 Monitors) + 1 MTRA Workstation + 1 Demo Workstation		
HIS/RIS Integration with Vendor:	Medistar		
Telecommunication with:	KH Wittlich, KH Bernkastel, KH Zell, KH Daun		
Special Applications or Integrations:	<p>Encrypted VPN over DSL (SDSL to 4 Hospitals, ADSL to 5 Homeoffices)</p> <p>PACS Server with Online-RAID 2,1TB – no Jukebox</p> <p>WEB-Distribution with a separate WEB-Server</p> <p>Expansion 2006 Interface CR, WEB in Bernkastel, Robot in Daun, Zell</p> <p>Server expansion</p>		

