



Hospital/Institution Name:	Donauries-Klinik		
City/Country:	Donauwörth	Germany	
Contact Persons (Phone, E-Mail):	Dr. Loitzsch	☎ +49 906 - 2999300	
Consultant (Address, Phone, E-Mail):	Med. Integral GmbH	☎ +49 8231 - 90832	
	Mr. Rupp	medhru@aol.com	
Connected Modalities DICOM:	1 x CT, 1 x MRI, 1 x Color Doubler, NUC		
Connected Modalities ANALOG:			
Installation date:	Feb. 2005		
Investment Volume ca.:	100'000 € - 120'000 €		
Patients per day:	GB per day:	45	1 GB
Server size:	TB per year:	1 TB	0,3 TB
Number of Diagnosis Stations:	6		
HIS/RIS Integration with Vendor:	RIS AKOSYSTEM interface		
Telecommunication with:	2 Hospitals		
Special Applications or Integrations:	Auto routing for picture to second hospital, full integrated workflow		





Donau Ries Hospital



Dr. Loitzsch Radiology Dr. Adolph Radiology Dr. Tscherne Nuc.M.



Turn able Diagnostic Workstation



Turn able Diagnostic Workstation