



Hospital/Institution Name: **GTK Polyclinic**

City/Country: Moscow Russia



Contact Persons (Phone, E-Mail): Dr. Samuriev Medical Director 


Consultant (Address, Phone, E-Mail): KBV Leipzig 


Connected Modalities DICOM: 1 x CT, 2 x CR, 3 x Laser Cam

Connected Modalities ANALOG: 1 x Mammo, 1 x Cathlab, 2 x US,

Installation date: 2005

Investment Volume ca.: € - €

Patients per day: GB per day: 210 GB

Server size: TB per year: 2,15 TB TB

Number of Diagnosis Stations: 6

HIS/RIS Integration with Vendor:

Telecommunication with: Telecommunication between two clinics and home offices

Special Applications or Integrations: GUI / User Manual etc. are provided in Russian Language.
 Old and new site of hospital complex are sharing Central Server

