


<b>Hospital/Institution Name:</b>		<b>Zollernalb Klinikum gGmbH, Albstadt and Balingen</b>	
<b>City/Country:</b>		Albstadt and Balingen	Germany 
<b>Contact Persons (Phone, E-Mail):</b>	Prof. Dr. med Michael Bitzer,	☎ +49 7431 99-1701	
	CA RAD	✉ michael.bitzer@zollernalb-klinikum.de	
<b>Contact Persons (Phone, E-Mail):</b>	Mr. Josef Weiss	☎ +49 7431 99-1000	
	GF	✉ josef.weiss@zollernalb-klinikum.de	
<b>Connected Modalities DICOM:</b>	2x CT, 1 x MRI (Praxis), 1 x CR, 3 x DR, 3 x US, 1 x Colour Doubler,		
<b>Connected Modalities ANALOG:</b>	1 x Video-Endoscopie (OP)		
<b>Installation date:</b>	April 2007		
<b>Investment Volume ca.:</b>	220 000,- € - € 240 000,-		
<b>Patients per day:</b>	<b>GB per day:</b>	5 GB	
<b>Server size:</b>	<b>TB per year:</b>	6 TB	1,5 TB
<b>Number of Diagnosis Stations:</b>	8		
<b>HIS/RIS Integration with Vendor:</b>	HIS and RIS (VIS) from iSoft,HL7-Interface to the VEPRO PACS VIS		
<b>Special Applications or Integrations:</b>	Clinic wide WEB image distribution through VEPRO's Image (WEB) Distribution Concept. Communication between the two Hospitals in the Hospital – group. VEPRO Information system (VIS), Worklist is coming from VEPRO VIS.		



Network Chart



Zollernalbklinikum Albstadt



Dr. med. Helmut Eier (OA RAD), Prof. Dr. med. Michael Bitzer (CA RAD)