


Hospital/Institution Name:	Diakoniekrankenhaus Halle		
City/Country:	Halle	Germany	
Contact Persons (Phone, E-Mail):	Herr Dr. Böhlke (Rad.)	+49 (0) 345 - 778 6 -0	
		+49 (0) 345 - 778 6 -399	
Consultant (Address, Phone, E-Mail):	Faust Consult GmbH Leipzig	+49 (0) 341 - 3936 -400	
	Herr Kast	leipzigfc@faust-consult.de	
Connected Modalities DICOM:	1 x CT, 2 x CT-WS, 1 x DSA, 1 x C-ARM, 2 x DLR, 1x MIC, 1x Laserimager		
Connected Modalities ANALOG:	1 x DL, 3 x Video		
Installation date:	December 2006		
Investment Volume ca.:	200.000 € - 250.000 €		
Patients per day:	GB per day:	60	3,7 GB
Server size:	TB per year:	2,4 TB	0,9 TB
Number of Diagnosis Stations:	7		
HIS/RIS Integration with Vendor:	PACS Interface to KIS (Boss)		
Telecommunication with:			
Special Applications or Integrations:	Clinic wide image distribution via DICOM from the start PACS server with online RAID 2,4 TB Redundant cluster server All modalities with DICOM worklist – interface Incl. 3 capsuled OP-WS with touch-screen Installation time: two weeks during rotate clinic-work-day		

References - Diakoniekrankenhaus Halle

