







Hospital/Institution Name:	DRK Krankenhaus Sömmerda		
City/Country:	Sömmerda	Germany	
Contact Persons (Phone, E-Mail):	Herr Streifler	+49 (0) 3632 - 67 - 0	
			✉ M.Streifler@khsdh.de
Consultant (Address, Phone, E-Mail):	Herr Krisch	+49 (0) 34671 - 65 163	
			✉ projekt1_mkbf@web.de
Connected Modalities DICOM:	1 x CT		
Connected Modalities ANALOG:	1 x Laserimager		
Installation date:	December 2006		
Investment Volume ca.:	50.000 € - 100.000 €		
Patients per day:	GB per day:	12	0,2 GB
Server size:	TB per year:	2,4 TB	51 GB
Number of Diagnosis Stations:	2		
HIS/RIS Integration with Vendor:	PACS Interface to KIS (medico//s) and RIS (Medavis)		
Telecommunication with:	3x Clinics (VPN-WAN 2 MBit symmetric)		
Special Applications or Integrations:	Permanent data distribution between 4 clinics – with VPN security Clinic wide image distribution via DICOM from the start PACS server with online RAID 2,4 TB Redundant cluster server All modalities with DICOM worklist – interface Installation time: one week during rotate clinic-work-day		

References - DRK Krankenhaus Sömmerda





DRK Krankenhaus Sömmerda