



Hospital/Institution Name:	Bezirkskrankenhaus Ebensfeld		
City/Country:	Ebensfeld	Germany	
Contact Persons (Phone, E-Mail):	Herr Drewitzki	+49 9547 81-2001	
			✉ Volker.Drewitzki@bezirksklinikum-obermain.de
Consultant (Address, Phone, E-Mail):	C-Consult	+49 731 69 203	
	Herr Christiansen		✉ christiansen@c-consulting.org
Connected Modalities DICOM:	CT, DR, Mobil DR, US, ES, C-Arm		
Connected Modalities ANALOG:	ES, C-Arm		
Installation date:	January 2008		
Investment Volume ca.:	480.000,00 € - 540.000,00 €		
Server size:	TB per year:	2,4 TB	0,6 TB
Number of Diagnosis Stations:	13		
Telecommunication with:	Image Center Dr. Bruckner Bamberg,		
Special Applications or Integrations:	Enterprise PACS Solution Interface to Siemens HIS (Medico//s) and RIS (WRAD). Interface to prosthetic planning from HecTec, all OT full integrated in the solution.		



Network Chart Hospital



Last edit date: 09.06.2010

Workflow Chart



Bezirksklinikum Obermain

