


Hospital/Institution Name:	Marienhospital Oelde / Sankt Elisabeth Hospital GmbH		
City/Country:	Oelde	Germany	
Contact Persons (Phone, E-Mail):	Herr Strake, IT-Manager	[49] 2522 99 4801	
			 l.strake@marienhospital-oelde.de
Consultant (Address, Phone, E-Mail):			
			
			
Connected Modalities DICOM:	Network Communication		
Connected Modalities ANALOG:			
Installation date:	March 2008		
Investment Volume ca.:	30.000,00 € - 80.000,00 €		
Patients per day:	GB per day:		
Server size:	TB per year:	2 TB	
Number of Diagnosis Stations:	3		
HIS/RIS Integration with Vendor:	Meierhofer AG. MCC		
Telecommunication with:	St. Elisabeth Hospital, Gütersloh		
Special Applications or Integrations:	VPN-Communication (DICOM send) Teleradiologie to St. Elisabeth KH Gütersloh		




Telecommunication Diagram



St. Elisabeth-Hospital Gütersloh



St. Elisabeth-Hospital Gütersloh



Testimonial Letter: