




Hospital/Institution Name:	Kreiskliniken Unterallgäu		
City/Country:	Ottobeuren, Mindelheim	Germany	
Contact Persons (Phone, E-Mail):	Herr Velter (Lt. IT)	+49 8261-797-7171	r.velter@kreisklinik-mindelheim.de
Consultant (Address, Phone, E-Mail):	Herr Trautner	+49 9571 7584320	e.j.trautner@sv-trautner.de
Connected Modalities DICOM:	CT, DSA, CR, C-Arm, Cathlab		
Connected Modalities ANALOG:			
Installation date:	August 2009 und December 2009 (Cardiologie)		
Investment Volume ca.:	180'000,00 € - 230'000,00 €		
Patients per day:	GB per day:	130	1,6
Server size:	TB per year:	1,2 TB	0,4 TB
Number of Diagnosis Stations:	5		
HIS/RIS Integration with Vendor:	Siemens MEDICO//S / Siemens WRAD		
Telecommunication with:	Klinikum Augsburg		
Special Applications or Integrations:	Satellite and main server solution. Full integrated interface to HIS and RIS per HL7 incl. reporting archiving in the PACS solution.		



Network Chart Hospital



HL7 technical diagram



Location Kreisklinikum Mindlheim



Location Kreisklinikum Ottobeuren