




Hospital/Institution Name:	Sankt Elisabeth Hospital GmbH		
City/Country:	Gütersloh	Germany	
Contact Persons (Phone, E-Mail):	Herr Fischer	+49 5241 507 5000	
			✉ Michael.Fischer@sankt-elisabeth-hospital.de
Consultant (Address, Phone, E-Mail):			
			✉
Connected Modalities DICOM:	CT, MRI, 3 x CR, Mammo, 2 x C-Arm, 6 x ES		
Connected Modalities ANALOG:			
Installation date:	06.06.2010		
Investment Volume ca.:	400'000 € - 450'000 €		
Patients per day:	GB per day:	250	3 GB
Server size:	TB per year:	15 TB	1,5 TB
Number of Diagnosis Stations:	12		
HIS/RIS Integration with Vendor:	Systema Compu Group		
Telecommunication with:	Image Center, more then 4 hospitals		
Special Applications or Integrations:	Cluster server solution. VIS for complete workflow management in the hospital. Full integrated interface to HIS. Interface to Praxis - system. Communication for Teleradiology (3 Hospitals) + Home office WS. Teleradiology workflow from hospital in Oelde.		




Network Chart Image centre



Interface



HL7 technical diagram



Workflow Mammography



Sankt Elisabeth Hospital