


|  |  |                       |   |
|--|--|-----------------------|---|
| <b>Hospital/Institution Name:</b>            | Radiologie Vorderpfalz   |                       |   |
| <b>City/Country:</b>                         | Ludwigshafen, Bad Dürkheim,<br>Frankenthal                                 | Germany               |  |
| <b>Contact Persons (Phone, E-Mail):</b>      | Dr. med. Uwe Himmighöfer   | +49 6322 94 17 79     |   |
|  |  | radiologe@t-online.de |   |
| <b>Consultant (Address, Phone, E-Mail):</b>  |  |                       |   |
| <b>Connected Modalities DICOM:</b>           | 1 x CT, 3 x MRI, NUK, Mammographie   |                       |   |
| <b>Connected Modalities ANALOG:</b>          |  |                       |   |
| <b>Installation date:</b>                    | First step 2001; Second step 2007, 3th step 2010                           |                       |   |
| <b>Investment Volume ca.:</b>                | 120.000 € - 180.000 €  |                       |   |
| <b>Patients per day: GB per day:</b>         | 110  |                       |   |
| <b>Server size: TB per year:</b>             | 5 TB   | 0,5 TB                |   |
| <b>Number of Diagnosis Stations:</b>         | 4  |                       |   |
| <b>HIS/RIS Integration with Vendor:</b>      |  |                       |   |
| <b>Telecommunication with:</b>               | Practice Ludwigshafen & 3 Home-Office-Workstations                         |                       |   |
| <b>Special Applications or Integrations:</b> | VPN- Connection to Practice Ludwigshafen and to 3 Home-office-workstations |                       |   |




**GESAMTLÖSUNG**  
Radiologie Ludwigshafen & Radiologie Bad Dürkheim & Frankenthal


**STANDORT – LUDWIGSHAFEN**



**STANDORT – BAD DÜRKHEIM**



**Neuanbindung: STANDORT – FRANKENTHAL**



Last edit date: 31.03.2010



Practice Entrance



DICOM – MRI