



Hospital/Institution Name:	Lukas Krankenhaus Bünde		
City/Country:	Germany		
Contact Persons (Phone, E-Mail):	Mr. Ulrich Falk	☎ +49 (0) 5223 / 167-0	
	www.lukas-krankenhaus.de	✉ u.falk@lukas-krankenhaus.de	
Consultant (Address, Phone, E-Mail):	Diagmed Consult, Mr. D. Riewe	☎ +49 (0) 2602 / 18626	
		✉ driewe@rz-online.de	
Connected Modalities DICOM:	MRI, Angio, 6x C-Arm, 2x CR, Lasercamera		
Connected Modalities ANALOG:	5 x X-Ray		
Installation date:	January 2011		
Investment Volume ca.:	400.000 – 500.000 €		
Patients per day:	GB per day:	200	5 GB
Server size:	TB per year:	16	1,5
Number of Diagnosis Stations:	13 + 13 Viewing Workstations		
HIS/RIS Integration with Vendor:	HIS: Siemens medico//s, RIS: Siemens WRAD		
Telecommunication with:	Paracelsus Hospital Gütersloh, 3 x home offices		
Special Applications or Integrations:	VIS back end, Worklist from VIS		
	10 x orthopaedic planning system		
	Cooperation with Lukas Hospital, Enger		

Network Chart



Last edit date: 04.05.2011

Workflow – Technical Diagram

