



Hospital/Institution Name:	Katharina-Kasper gGmbH „St. Marienkrankenhaus“ & „St. Elisabeth Krankenhaus“		
City/Country:	Frankfurt	Germany	
Contact Persons (Phone, E-Mail):	Mr. Walter Semmler, CEO	+49 (0) 69 / 1563-0	
	http://www.katharina-kasper.de/	w.semmler@katharina-kasper.de	
Consultant (Address, Phone, E-Mail):	Mr. Hagen Woecht	+49 (0) 30 / 293618-0	
	Ingenieurgesellschaft Mengel & Partner GbR	woecht@mengelundpartner.de	
Connected Modalities DICOM:	St. Marien KH: CT, C-Arm, CR-Systems St. Elisabeth KH: , CT, Uroscopy, CR-Systems		
Connected Modalities ANALOG:	St. Marien KH: Fluoro, 3x Ray, US, 5 x Endo, C-Arm, Gamma-Camera, St. Elisabeth KH: Fluoroscopy, 2 x C-Arm, 5 x Endoscopy		
Installation date:	April 2011		
Investment Volume ca.:	1.000.000 € - 1.200.000 €		
Patients per day:	GB per day:	225	4 -5
Server size:	TB per year:	25 TB	1 – 1,5 TB
Number of Diagnosis Stations:	St. Marien KH: 16; St. Elisabeth KH: 7		
HIS/RIS Integration with Vendor:	KIS: Tieto iMedOne		
Telecommunication with:	Radiologic practise Adickesallee, MVZ St. Kamilus GmbH,		
Special Applications or Integrations:	VIS in Frontend incl. Ordering, Schedule, Documentation, X-Ray Documentation, Reporting, Statistics Oncologic treatment system & Oncologic planning system Integration of the radiology department, operating room, nuclear medicine, radiation therapy		


Network Chart St. Marien KH




Network – St. Elisabeth KH



Interfacing – Technical chart



VEPRO Katharina-Kasper gGmbH
RIS / PACS Workflowkonzept



Last edit date: 11.03.2010

Workflow chart



St. Marienkrankenhaus Frankfurt/Main



St. Elisabethen-Krankenhaus Frankfurt/Main