








Hospital/Institution Name: Hospitaux Gabonais, Agondje,			
City/Country:		Gabon	Libreville 
Contact Persons (Phone, E-Mail):		M.SC. Eng. Martin Schönthaler 	
Consultant (Address, Phone, E-Mail):		VAMED ENGINEERING GmbH & CO KG	 +43-1-601727 227
		M.SC. Eng. Martin Schönthaler	 martin.schoenthaler@vamed.com
Connected Modalities DICOM:		2 x CT, Orthopantomograph, NUC, Fluoroscopy, 2 x C-Arm, 4 x Ultrasound, 2 x Lasercamera; Radiotherapy-Treatment Planning System, 2 x CR	
Connected Modalities ANALOG:		4 x Endoscopy	
Installation date:		November 2011	
Investment Volume ca.:		350.000 €- 400.000 €	
Patients per day:	GB per day:	400 P / day expected	5 GB
Server size:	TB per year:	5 TB	1,4 TB
Number of Diagnosis Stations:		2	
HIS/RIS Integration with Vendor:		-	
Telecommunication with:		-	
Special Applications or Integrations:		VEPRO Information System (VIS) is the leading system and the workflow is fully operated by VEPRO systems. VIS organized, in- and outpatient registration, ordering, scheduling, reporting, documentation and controls the VEPRO PACS. Also have the hospital with the VEPRO System an electronic medical record.	

Network Chart



Last edit date: 20.01.2012

Network – Technical Diagram





The Hospital



Entrance Hall




Inpatient Buildings



Block of buildings



Hospital area



Training reading Physicians