




<b>Hospital/Institution Name:</b>	Asklepios Klinik Melsungen		
<b>City/Country:</b>	Melsungen	Germany	
<b>Contact Persons (Phone, E-Mail):</b>	Hr. Exnowski	+49 56 61 7 77 0	
			✉ m.exnowski@asklepios.com
<b>Consultant (Address, Phone, E-Mail):</b>			
			✉
<b>Connected Modalities DICOM:</b>	1 x CT, 1 x CR, 1 x MWS		
<b>Connected Modalities ANALOG:</b>			
<b>Installation date:</b>	January 2011		
<b>Investment Volume ca.:</b>	100.000,00 € - 130.000,00 €		
<b>Server size:</b>	<b>TB per year:</b>	4,0 TB	0,5 TB
<b>Number of Diagnosis Stations:</b>	2		
<b>HIS/RIS Integration with Vendor:</b>	AGFA Orbis		
<b>Telecommunication with:</b>			
<b>Special Applications or Integrations:</b>	The solution is part of a satellite server concept from 3 locations (first in Schwalmstadt, second in Bad Wildungen, third in Melsungen)		



Network Chart

**VEPRO** Asklepios Klinikum Schwalmstadt Standort Melsungen



Last edit date: 31.03.2011

HL7 – Technical Diagram



Asklepios Klinik Melsungen