

# VEPRO eHEALTH Solutions



WEBstudio Server



## WEBstudio Server – for small and midrange Hospitals

One Hardware – One installation – Any own device

- Full PACS or RIS – scalable Software to customer needs
- Scalable from 10 to 50 users
- Connection of any existing own devices: Desktop or Mobile
- Access inhouse Hospital via LAN or WLAN or from any location by Internet
- Easy installation and full operation in two days
- Incl. standard service 8:00-18:00 – incl. Updates  
Standard Warranty: 3 years

### Technical specifications:

- Intel® Xeon® - 6Core CPU
- 64GB RAM
- 2x RJ45, 1000MBit Network on board (CU)
- 2x Harddrive 4000 GB RAID1 or 4x Hard drive, 4000GB, SATA, RAID6
- Operating System:
  - Debian 9 with Proxmox Virtual Environment 5.4
  - Windows® Server 2019 Standar