


## Perfect for your reception:

- ❖ Compact
- ❖ Smart
- ❖ Reliable



Just let DICOM CDs import by the patients themselves...  
insert CD - wait - remove CD!

### Features

- Import all DICOM Patient CDs / DVDs
- Storage of patient data in a local FIFO archive
- Send the data directly to any PACS
- Display of patient studies via the integrated DICOM viewer

### Options

- (G) DICOM Worklist SCU - EMR License - VEPRO PACS Broker

### Technical specifications

<b>Mini PC:</b>	CPU – Intel Celeron J4005, 4GB DDR RAM, 2x HDMI, 250GB SSD Network: 1 GB RJ45, Bluetooth, WLAN 4 x USB 3.0, 2 x USB 2.0 Memory Card slot, Audio, 110/220 V – 115x111x51 mm Windows® 10 IoT Ent 2019 LTSC, 64 Bit
<b>CD / DVD drive:</b>	Slot-In, Speed 10x – up to 5Gbit/s – USB 2.0 - 142x138x16 mm
<b>Monitor:</b>	10,1 inch/25cm – Touchscreen – 1280x800 – HDMI, VGA – USB 2.0 – Audio 110/220 V – 253x166x35 mm
<b>Software:</b>	VEPRO EMR-Manager (PACS) CE 0297 - DIN EN ISO 13485 : 2016 + AC : 2017-07 Base Viewer <a href="https://www.vepro.com/en/pacs.html">https://www.vepro.com/en/pacs.html</a>  (U) EMR License, DICOM Send Servers - SCU Add on module to send DICOM data to any 3rd party DICOM Node Languages: German, English, French, Spanish, Russian